


# Series 8012, 8020, 8040, 8050 and 8060 Indicating Moisture Traps



These units are designed to remove water, oil and organics from gases commonly used as gas chromatography carrier gases. They are constructed from Lexan® polycarbonate tubing with aluminum end caps sealed with Viton® o-rings, except for the 8050 which has a solid aluminum housing and is thus non-indicating. All units are filled with a mixture of molecular sieve 13X and indicating molecular sieve 4A. These are the highest capacity molecular sieves available and the preferred choice for gas drying. The blue indicating sieves turn buff color at 20% relative humidity.

## Features

- Reduces Water to Less than 20 ppb
- Inlet and Outlet O-ring Sealed  
Connectors are Equipped with 40 Micron Sintered Stainless Steel Frits to Prevent Particulates from Entering Your System (Model 8060 is 100 Micron)
- Mounting Clip Available for Convenient Installation

## Specifications

Moisture Removal Capacity . . . . .	
8012	21.6 grams
8020	36.0 grams
8040	72.0 grams
8050	132 grams
8060	245 grams
Dimensions . . . . .	
8012	1.5" O.D. x 9.0" long
8020	1.5" O.D. x 12.5" long
8040	1.75" O.D. x 17.5" long
8050	2 3/8" O.D. x 17" long
8060	3" O.D. x 19" long
Working Pressure . . . . .	
8012, 8020, 8040, 8060	125 PSIG
8050	250 PSIG
High Flow Capacity . . . . .	
8060	up to 20 SCFM

Gas traps should be mounted in a vertical position to ensure proper contact of the gas with the adsorbent. Use model 8012C, 8040C, 8050C, or 8060C mounting clip with 8012, 8020, 8040, 8050 and 8060 Series moisture trap.

## How to Order:

Model Number	Capacity	Connections
8012-2*	120 cc	1/8" tubing compression fittings
8012-4*	120 cc	1/4" tubing compression fittings
8020-2*	200 cc	1/8" tubing compression fittings
8020-4*	200 cc	1/4" tubing compression fittings
8040-2*	400 cc	1/8" tubing compression fittings
8040-4*	400 cc	1/4" tubing compression fittings
8040R	400 cc	Provides enough for three 120 cc, two 200 cc, or one 400 cc refill
8050-8*+	735 cc	1/2" tubing compression fitting
8050R	1500 cc	Provides enough for two refill
8012C		Clip for mounting 8012 and 8020 units
8040C		Clip for mounting 8040 units only
8050C		Clip for mounting 8050 units only
8060-8SS	1350 cc	1/2" tubing stainless steel compression fittings
8060R		Provides enough molecular sieve mix for two refills
8060C		Mounting clip (two per unit suggested)

+8050 is a non-indicating trap

\*Available with stainless steel compression fittings - add "SS" to part number