


Liquid Nitrogen Handling and Safety Accessories



Stainless Steel Transfer Dewars

Used for working safely with small amounts of liquid nitrogen. NEVER use thermos bottles, Styrofoam coolers, or other unapproved containers.

Size	Part No.
5 Liter	11014070
2.5 Liter	11048940
1.5 Liter	11053301



Cryo Apron

Part No.
10464394



Cryo Gloves


Size	Part No.
Mid-arm Length, Medium	9717119
Mid-arm Length, Large	9717129
Mid-arm Length, Extra Large	9717139
Elbow Length, Medium	9717149
Elbow Length, Large	9717159
Elbow Length, Extra Large	9717169



More details on page 228

Face Shields

Style	Part No.
Heavy Duty Dual Crown Face Protection	61001400
General Purpose Face Protection	Crown: 61001440 Window: 61001300



Cryogenic Transfer Hose

1/2" ID Stainless steel, flexible hose for safely transferring liquid nitrogen.

Size	Part No.
4' Length	9713159
6' Length	9713109



Phase Separator

Reduces flash loss and hazardous splashing during liquid nitrogen transfer to storage dewars.

Size	Part No.
Small, 11/16" diameter x 1 1/2" long	10615885
Large, 1 5/16" diameter x 3" long	10615877
Extra Large, 1 5/8" diameter x 4" long	10615869



T-Valve Assembly

Allows two Cryogenic Transfer Hoses to be safely attached to one liquid nitrogen supply source.

Part No.
10784443